# RegenOss Ortho € 0373

---

Italian IT

English EN

French FR

German DE

Spanish ES

Portuguese PT

Turkish TR

Dutch NL

Bulgarian BG

Slovak SK

Czech CZ

Slovenian SLO

---

	<p>Fabbricante  Manufacturer  Fabricant  Hersteller  Fabricante  Fabricante  Üretici  Fabrikant  Производител  Výrobca  Výrobce  Proizvajalec</p>
	<p>Non risterilizzare  Do not resterilize  Ne pas restériliser  Nicht reesterilisieren  No volver a esterilizar  Não submeter à reesterilização  Sterilize etmeyiniz  Niet opnieuw steriliseren  Да не се стерилизира повторно  Nesterilizujte opakovane  Ne re-sterilizovat  Ne sterilizirajte znova</p>
	<p>Non riutilizzare  Do not re-use  Ne pas réutiliser  Nicht wiederverwenden  No volver a usar  Não reutilizar  Yeniden kullanmayınız  Niet hergebruiken  Да не се използва повторно  Neroužívejte opakovane  Pouze na jednorázové užití  Samo za enkratno uporabo</p>
	<p>Non utilizzare se la confezione è danneggiata  Do not use if the package is damaged  Ne pas utiliser si l'emballage est endommagé  Nicht benutzen, wenn die Verpackung beschädigt ist  No se debe utilizar si el envase está dañado  Não utilizar se a embalagem estiver danificada  Ambalaj hasarlı ise kullanmayınız.  Niet gebruiken als de verpakking beschadigd is  Да не се употребява, ако опаковката е повредена!  Neroužívejte, ak je obal poškozený  Při poškozeném obalu nelze použít  Izdelka ne uporabljajte, če je embalaža poškodovana</p>

	<p>Sterile, incluso il metodo di sterilizzazione (radiazioni ionizzanti)  Sterile, including the method of sterilization (ionizing radiations)  Stérile, incluant la méthode de stérilisation (radiations ionisantes)  Steril, samt Sterilisationsmethode (ionisierende Strahlungen)  Estéril, incluido el método de esterilización (radiaciones ionizantes)  Estéril, incluindo o método de esterilização (radiações ionizantes)  Steril, sterilizasyon yöntemi (iyonlaştırıcı radyasyon) dahil  Steriel, inclusief sterilisatiemethode (ioniserende straling)  Стерилен, включително метод на стерилизация (ионизиращо лъчение)  Sterilné, metóda sterilizácie (ionizujúce žiarenie)  Sterilní, včetně sterilizační metody (ionizující záření)  Sterilno, vključno z načinom sterilizacije (ionizirajoče sevanje)</p>
	<p>Utilizzare entro  Use by  Date limite d'utilisation  Mindestens haltbar bis  Fecha de caducidad  Usar antes de  Son kullanma tarihi  Te gebruiken tot  Используй до  Dátum expirácie  Datum použitelnosti  Rok uporabnosti</p>
	<p>Marcatura CE  CE marking  Marquage CE  CE-Kennzeichnung  Marcado CE  Marcação CE  CE işareti  CE-markering  CE маркировка  CE značka notifikovaného organu  CE značka  Oznaka CE</p>
	<p>Numero di lotto  Batch code  Numéro de lot  Chargennummer  Número de lote  Número de lote  Parti numarası  Lotnummer  Партиден №  Číslo šarže  Číslo šarže  Številka serije</p>

	<p>Numero di catalogo          Catalogue number          Numéro de catalogue          Katalog-Nummer          Número catálogo          Número de catálogo          Katalog numarası          Catalogusnummer          Каталоген №          Katalógové číslo          Katalogové číslo          Kataloška številka</p>
	<p>Numero di pezzi          Number of pieces          Nombre de pièces          Stückzahl          Número de piezas          Número de unidades          Adet sayısı          Aantal stuks          Бройки от продукта          Počet kusov v balení          Počet kusů          Število kosov</p>
	<p>Attenzione, vedere le istruzioni per l'uso          Warning: see instructions for use          Attention, voir les instructions pour l'emploi          Achtung, siehe die Gebrauchshinweise          Cuidado, consultar las instrucciones de uso          Atenção, consultar as instruções de uso          Dikkat, kullanım talimatlarına bakınız          Let op, zie de gebruiksaanwijzing          Внимание, прочетете инструкцията за употреба!          Pozor, přeštudujte si příložený návod na použití          Pozor, přečtěte si návod k použití          Opozorilo: upoštevajte navodila za uporabo</p>
	<p>Consultare le istruzioni operative          Consult the operating instructions          Voir les instructions pour l'utilisation          Bedienungsanleitung beachten          Consultar las instrucciones operativas          Consultar as instruções operacionais          İşletim talimatlarına başvurunuz          Lees de gebruiksaanwijzing          Консултирайте се с инструкцията за употреба!          Dodržujte návod na použití          Prostudujte si pokyny          Preverite navodila za uporabo</p>

# RegenOss - Ortho

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### DESCRIZIONE

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## Bone substitutes

### DESCRIPTION

RegenOss is a medical device for bone defect reconstruction.

RegenOss is an implantable, reabsorbable, non-active medical device, acting as bone substitute. It is composed of type-I collagen (of equine origin) and biomimetic hydroxyapatite whose chemical, chemical-physical and biological features are equivalent to those of human bone mineral constituent.

Collagen is a fibrous protein and represents the most important structural component of many human connective tissue extracellular matrixes, while hydroxyapatite is the main bone mineral element (65% by weight).

The presence of type-I collagen makes this device chemically and structurally biomimetic and ensures great in-site stability.

RegenOss is a flexible composite that can reproduce the bone anatomical structure. RegenOss, thanks to its special configuration and biomimetic chemical composition, can promote bone tissue regeneration processes.

### INDICATIONS for USE

RegenOss is indicated for recovering bone substance loss, in particular:

- in presence of bone defects due to trauma, degenerative diseases of the skeletal system or surgical excisions, in combination with osteosynthesis implants and/or prosthesis;
- for filling defects due to autologous bone removal or to tumour and cyst removals;
- for treating pseudarthrosis;
- in presence of bone deficiencies following prosthetic implant revisions;
- in all pathological conditions in the orthopaedic surgery field that require the generation of new bone tissue.

RegenOss can also be used mixed with antibiotics.

The device acts as a bone substitute that is reabsorbed during the remodelling process of the newly-formed tissue.

RegenOss can be used in combination with saline fluids to improve its flexibility. The device can be used with surgical glues to increase stability in case of application on tissue defects with uneven or jagged surfaces.

### AVAILABLE SIZES

The device is available in various pre-shaped forms and sizes.

## INSTRUCTIONS for USE

Instructions for the surgical use of RegenOss are described below.

The following instructions are to be considered only as recommendations for the surgical use of RegenOss, they are not intended to modify or replace standard procedures for the treatment of bone defects. Surgical procedures involving the use of bone substitutes may have variable outcomes. The factors to be taken into account when choosing a bone substitute and the surgical technique are as follows: age and bone quality of the patient, location of the defect, expected loading condition, level and quality of tissue vascularity near the defect, possibility for the substitute to be in direct contact with the vital bone tissue adjacent to the defect, chances to ensure the complete filling of the bone defect, possibility to obtain correct and sufficient primary stabilisation of the bone substitute, chances to obtain complete coverage of the bone substitute.

### **Pre-operative treatment**

It is advisable to carry out pre-operative radiographs to evaluate the conditions of the bone segment and the size of the defect to be treated as well as to choose the type and amount of product to apply and the fixation methods and devices. In the event the surgeon decides to proceed with surgery, the patient must undergo antibiotic prophylaxis. An alternative prophylaxis is essential to be used in patients allergic to penicillins.

It is important to carefully verify the absence of infections or inflammations prior to surgery.

### **Intra-operative aspects**

The device must be used only for surgeries performed in sterile operating environments, adopting every precaution. The package is to be opened after verifying its integrity and only in sterile environments.

RegenOss does not interfere with normal surgical procedures. After careful cleaning of the bone surface, the product can be directly applied onto the bone area to be treated.

If necessary, shape the device according to the site by means of sterile forceps or scalpels. The quantity of device to be used must not exceed the size of the defect.

Place the device in the bone defect to be treated and stabilise it, pressing gently.

If necessary, use surgical glue for further stabilisation.

It is essential to prevent any dispersion of material outside the implant site but the amount of product used must be enough to fill the entire bone defect.

Avoid irrigating the defect site once the substitute has been positioned.

Excess material must be carefully removed making sure not to dislocate the bone substitute placed in the defect.

Once the filling is completed, close the surgical site following the standard procedure.

Bleeding of the bone compartment may facilitate stabilisation of the device.

Close the surgical site according to standard procedures.

RegenOss does not have sufficient mechanical resistance to bear direct loads in the immediate post-op phase. If applied in areas that require immediate loading, it is important to use suitable fixation means or synthetic systems.

### **Post-operative treatment**

The patient should undergo a suitable post-operative antibiotic therapy according to the surgical procedures.

The surgeon must give the patient all the instructions related to the specific application to guarantee a proper post-operative course.

In absence of complications and if all the surgeon's post-op indications are followed correctly, primary stabilisation of the implant generally occurs within 6-8 weeks thanks to the osteointegration process. However, it is important to verify through radiographs the newly formed bone tissue before subjecting the area to direct loading.

## **CONTRAINDICATIONS**

Do not use the device in case of:

- systemic or localised infection;
- inflammatory or autoimmune disease;
- hypercalcemia;
- degenerative bone disease;
- coagulation disorders;
- metabolic disorders;
- diabetes;
- thyroid function abnormalities or complications;
- systemic conditions affecting wound healing;
- allergies to collagen of equine origin and to calcium phosphate salts.

The product must not be implanted in areas where complete coverage with soft tissue is not possible.

## **USE IN SPECIFIC SUBJECTS**

RegenOss efficiency is unknown in case of:

- renal diseases;
- pregnancy;

- breast-feeding;
- ongoing radiotherapy;
- neoplasia;
- cardiovascular diseases.

No side effects due to the use RegenOss in presence of cancer have been reported. However, patients undergoing bone filler treatment following tumour removal must be carefully monitored.

## WARNINGS

RegenOss is intended to be used only by qualified and specialised medical personnel.

RegenOss does not have sufficient mechanical resistance to bear direct loads before regeneration of the newly formed tissue. If applied in areas that require immediate stabilisation and/or loading, it is important to use fixation means or synthetic systems.

Do not use the product after its expiry date.

RegenOss is a single-use product. Do not use any remains for other operations or other patients but dispose of it as contaminated hospital waste according to the regulations in force at the hospital facility. The product must not be re-sterilised.

The manufacturer is not liable for risks deriving from processes to recondition the device; reconditioning of the device implies the loss of its features and functions.

The product is not to be used if its packaging is damaged. The sterility of the product is guaranteed only if the package is unopened and undamaged.

The package contains an additional label displaying product tracking information; the implant must be tracked with the label, which must be applied to the patient's medical records.

In case of allergic reaction, the device must be removed and the patient must be treated with an adequate therapy.

The factors to be taken into account when choosing the device and the surgical technique are as follows:

- age of the patient;
- bone tissue quality;
- localisation of the defect;
- loading conditions.

Surgical procedures involving the use of bone substitutes may have variable outcomes, which also depend on the patient's compliance with the post-op protocol.

## PRECAUTIONS

In order not to compromise the porosity of the product, do not use an excessive amount of material or press the device too much when placing it in the defect site. RegenOss does not guarantee immediate post-implantation mechanical support therefore, it is important to use an adequate fixation system and walking aid.

## STERILIZATION

All RegenOss formats are supplied sterile (sterilised by 25 kGy gamma-rays). The device is a single-use product and cannot be resterilised. Its reuse, in whole or in part, may involve the risk of cross-contamination and the risk of infection at the implant site.

## STORAGE

Keep the product in a cool and dry place, protected from light and from heat sources (+10 °C / +30 °C).

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**FIN-CERAMICA FAENZA S.p.A.**