**JRI ORTHOPAEDICS LTD**  
18 Churchill Way,  
31A Business Park,  
Chapelton,  
Sheffield,  
S35 2PY, UK

Instructions for Use

**JRI 155-023**

**Furlong<sup>®</sup> Modular Total Hip  
Replacement System**

# Furlong® Modular Total Hip Replacement System

## Contents

---

English .....	3
Deutsch .....	XX
Français .....	XX
Italiano .....	XX
Español.....	XX
Português .....	XX
Norsk.....	XX
Svenska .....	XX
Suomi .....	XX
Dansk .....	XX
Nederlands .....	XX
Polski.....	XX
Česky .....	XX

## Important Information

Please read instructions for use and corresponding operative technique, prior to use in a clinical setting. The Surgeon should be familiar with the appropriate operative technique.

## Material(s)

Material(s)		Standards
High Nitrogen Stainless Steel	HNSS	ISO 5832-9
Cobalt Chrome	CoCr	ISO 5832-4
Ultra-High Molecular Weight Polyethylene	UHMWPE	ISO 5834-2
Cross-linked Ultra-High Molecular Weight Polyethylene	XL-UHMWPE	ISO 5834
BIOLOX® <i>delta</i> ceramic		ISO 6474-2

## Symbols

Heads: S = Short neck M = Medium neck L = Long neck XL = Extra Long neck

Stems: LS = Long Spigot 5 = 5mm lateral offset XXS = Extra Extra Small. XS = Extra Small.  
S = Small. M = Medium L = Large

## Intended Use

The **Furlong® Modular Total Hip Replacement System** is a multi-component system designed for use in cemented primary or revision total hip arthroplasty

## Treatment Population

**Furlong® Modular Total Hip Replacement System** is designed for use in skeletally mature individuals with sufficient bonestock to support implantation and fixation of the device.

## Description

The following components comprise the **Furlong® Modular Total Hip Replacement System**:

The **Furlong® Modular Femoral Stem**: a highly-polished, straight double tapered, collarless prosthesis, made of high nitrogen stainless steel and for use with cement. These are available in a variety of sizes (XXS, XS, S, M, and L); in both 137° and 140° neck angles; and in both standard and 5 mm lateral offset.

The **Furlong® Modular Femoral Revision Stem**: a highly-polished, straight double tapered, collarless prosthesis, made of high nitrogen stainless steel and for use with cement (same design as the standard stem). These have a 140° neck angle, and are available in a variety of sizes (S, M, and L); and in both standard and 5 mm lateral offset. Two stem lengths are available: 200 mm and 300 mm; the longer distal length is designed to provide additional stability. These stems are for use for revision total hip replacement only.

**CoCr & BIOLOX®*delta* heads**: modular femoral heads with a 12/14 taper; and various materials, sizes and neck lengths are available for use with this system. **Cobalt Chrome heads** can be used with polyethylene acetabular liners only; **BIOLOX®*delta* ceramic heads** can be used with both polyethylene and ceramic liners. Special Instructions for Use are required for ceramic heads: Refer to JRI Instructions for Use 155-020

The **JRI Muller Type Acetabular Cup**: a cemented Muller type acetabular cup available in a variety of sizes for use with cement. They are made of UHMWPE with an integral radiopaque stainless steel X-ray wire; the larger bore sizes are available in highly cross-linked UHMWPE.

**Femoral Cement Plugs**: made of UHMWPE with an integral radiopaque stainless steel X-ray wire and available in a variety of sizes. They may be used to contain the pressurised cement in the femoral canal at the correct depth.

# Furlong® Modular Total Hip Replacement System

English

---

Extensive clinical use of this design since 1979 has proven its biomechanical stability and biocompatibility. Various sizes of implants are available to accommodate anatomical variations. Smaller sized implants are intended for patients with small bone and normally slight weight and could be inappropriate for other patients.

## Note

The **Furlong® Modular Total Hip Replacement System** can be used with other compatible components of the JRI Furlong® System only, with the corresponding taper connection.

## Indications

The **Furlong® Modular Total Hip Replacement System** is indicated for use in the following conditions:

- Severely painful and/or disabled joint from osteoarthritis, traumatic arthritis, rheumatoid arthritis, or congenital hip dysplasia
- Avascular necrosis of the femoral head
- Acute traumatic fracture of the femoral head or neck
- Failed previous hip surgery including joint reconstruction, internal fixation, arthrodesis, hemi-arthroplasty, surface replacement arthroplasty, or total hip replacement
- Certain cases of ankylosis

## Contra-indications

The device should NOT be implanted where there is active infection, insufficient bonestock to support the prosthesis or provide adequate fixation.

Further contra-indications may include the following conditions:

- Severe muscle, nerve or vascular diseases that endanger the extremity in question
- Severe deformations, tumours
- Severe osteoporosis or deficient bone substance that may endanger stable seating of the prosthesis
- Systematic and metabolic disorders
- Weakened or compromised immune system (HIV, tumours, infections)

## Precautions

The following conditions require caution and due consideration during pre-operative planning by the surgeon:

- Obese or severely overweight patients
- Excessive loading through arduous activity
- Lack of mental faculties to understand post-operative recuperative regime
- Alcohol dependency or drug abuse
- History of falls or disabilities

In patients' with a high body mass index, where delayed surgery is feasible, it is advisable that a programme of weight reduction is undertaken prior to any joint replacement surgery.

## Warning

**This device should only be implanted by operating surgeons who are familiar with the general problems of prosthetic surgery and who are sufficiently trained to the product-specific operative technique. The surgeon is responsible for ensuring the surgery is carried out properly and in line with the instructions provided in the operative technique.**

**As a manufacturer, JRI Orthopaedics Ltd. Is not responsible for any complications arising from incorrect diagnosis, incorrect choice of implant, incorrect operative technique, treatment methods limitations or inadequate asepsis.**

## Pre-operative

Pre-operative planning allows the surgeon to assess for implant size and restoration of biomechanics. X-ray templates should be used to estimate implant sizes, placement and joint alignment. Failure to carry out proper planning may lead to incorrect choice of implant type/size. Ensure all implant sizes and required instrumentation are available prior to surgery. Consult operative technique and training materials provided by JRI Orthopaedics Ltd. before use.

Functionality of surgical instruments should be checked. Use of damaged instruments may lead to early failure of implant.

Transient bacteraemia can occur after surgical procedures. To prevent late infection at the implant site, many orthopaedic surgeons advise the use of antibiotic prophylaxis before and after such procedures for their patients.

## Intra-Operative

Always use a trial for any test fit and to check the range of motion. Failure to use the optimum size of implant of ensuring sufficient fixation may result in early revision surgery.

Care should be taken not to over tighten bone screws. Ensure the appropriate selection of bone screw length and location to avoid damage to underlying soft tissue areas.

Implants **MUST NOT** be re-used under any circumstances as the functionality, integrity and/or sterility of that device may have been adversely affected and therefore cannot be guaranteed. Implant components from JRI Orthopaedics Ltd. should not be used together with those of another manufacturer, since compatibility of mating parts cannot be assured.

## Post-Operative

Physicians should ensure patients are aware of implant loading limitations and ensure consistent post-operative care, from a suitably qualified professional, is made available. The incidence and severity of complications are usually greater in surgical revisions than primary operations. Early detection of an impending complication allows for timely and effective countermeasures.

Post-operative care should incorporate recognised procedures and take into account information from the operative technique and documented according to internal hospital procedures. Failure to do this may result in malalignment, delayed wound/bone healing, implant failure, infection or impaired joint function.

## Side Effects

As with all major surgical procedures, side effects and adverse events can occur. Some of the more common complications include:

- Problems resulting from anaesthesia and patient positioning (pain, nausea etc.)
- Problems resulting from use of bone cement (i.e. Bone Cement Implantation Syndrome)
- Infection
- Aseptic loosening of implant
- Dislocation, subluxation, insufficient range of movement, leg length discrepancy
- Damage to soft tissue
- Fracture of implant, bone or cement
- Allergy/hypersensitivity reactions
- Cardiovascular/pulmonary embolism, respiratory infection, venous thrombosis, neuronal dysfunction, haematoma or delayed wound healing

# Furlong® Modular Total Hip Replacement System

English

## Revision

A range of **Furlong® Modular Femoral Revision Stem** based on this system, are available and larger sizes are offered for revision applications.

Ensure all fragments of the primary prosthesis and any bone cement (if applicable) are removed. Clean and prepare the area for implantation in line with current Instructions for Use and Operative Technique of selected revision implant, prior to insertion of device.

Special instructions apply when ceramic components are revised. If revision of any ceramic component, always use a ceramic formal head. If revision of a ceramic femoral head and the femoral stem remains in-situ, always replace with a dedicated JRI revision ceramic head with a Ti taper sleeve, BIOLOX® OPTION. **Never** use another ceramic head without a taper sleeve on an existing femoral stem or damaged taper. **Never** use a metal femoral head following the revision of any ceramic component. If revision is required in case of a fracture of a ceramic component, make sure all ceramic particles have been removed. Refer to JRI Instructions for Use 155-020 (ceramic components) prior to revision surgery for more details.

## Storage & Handling

Store implants in their original protective packaging in a clean and dry atmosphere.

Do not use if the packaging is damaged. Do not use this product after the expiry date (year-month) shown on the product packaging.

Avoid removing from packaging until immediately prior to use. Inspect device prior to use. Visibly damaged, scratched, improperly handled implants and implants that have already been used must not be implanted under any circumstances as the functionality, integrity and/or sterility of that device may have been adversely affected and therefore cannot be guaranteed. Exposed articular surfaces must be neither marked nor come into contact with metallic or hard objects. Touching the articular surfaces, the Hydroxyapatite coating or the tapered interface for the liner must be avoided.

## MRI Safety Information

Non-clinical testing has demonstrated that **Furlong® Modular Total Hip Replacement System** is MRI Conditional. A patient with this device can be safely scanned in an MR system meeting the following conditions:

- Static magnetic field of 3.0 T
- Maximum spatial field gradient of 720 gauss/cm
- Maximum whole body averaged specific absorption rate (SAR) of 3.0W/kg

Under the scan conditions defined above, the Furlong® Modular THR are expected to produce a maximum temperature rise of 1.7°C after 15 minutes of continuous scanning.

The image quality may be compromised if the area of interest is in the exact same area, or relatively close to, the position of the medical device. Therefore, optimization of MR imaging parameters to compensate for the presence of this device may be necessary:

Pulse Sequence	T1-SE	T1-SE	GRE	GRE
Signal Void Size	42,763-mm <sup>2</sup>	17,875-mm <sup>2</sup>	56,249-mm <sup>2</sup>	28,917-mm <sup>2</sup>
Plane Orientation	Parallel	Perpendicular	Parallel	Perpendicular

Additional information regarding MRI safety is available upon request at JRI Orthopaedics Ltd.

# Furlong<sup>®</sup> Modular Total Hip Replacement System

English

---

## **Functional Device Lifespan**

The functional lifespan of an implant may be impacted by surgeon and selected operative technique, patient physiology and activity levels. Whilst it is normally expected that the lifetime of this implant will exceed a minimum of 10 years, it will be subject to wear and tear through normal everyday use.

## **Products intended for single-use must not be re-used**

Re-use or reprocessing (e.g. cleaning and re-sterilisation) may compromise the structural integrity of the device and / or lead to device failure which may result in patient injury or death.

Reuse or reprocessing of single-use devices may create a risk of contamination (e.g. transmission of infectious material) which may result in injury or death.

## **Further information**

For product specific operative technique training or further information, please contact your JRI Orthopaedics Ltd. Sales Representative or JRI Orthopaedics Ltd. directly.

**JRI ORTHOPAEDICS LTD**

18 Churchill Way,  
35A Business Park,  
Chapelton,  
Sheffield,  
S35 2PY, UK

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Další informace o bezpečnosti MRI vám na požádání poskytne společnost JRI Orthopaedics Ltd.