JRI ORTHOPAEDICS LTD
18 Churchill Way,
31A Business Park,
Chapelton,
Sheffield,
S35 2PY, UK

Instructions for Use

JRI ACE Acetabular Cup System

English	3
Czech	9
Dansk	15
Deutsch	21
Español	28
Suomi	35
Français	41
Italiano	48
Netherlands	55
Norsk	62
Polski	68
Português	75
Svenska	82

Important Information

Please read instructions for use and corresponding operative technique, prior to use in a clinical setting. The Surgeon should be familiar with the appropriate operative technique.

Material(s)

Material(s)		Standards
Titanium Alloy	Ti6Al4V	ISO 5832-3
Hydroxyapatite	CA ₅ OH(PO ₄) ₃	ISO 13779-1
Cobalt Chrome	CoCr	ISO 5832-4
Cross-linked Ultra-High Molecular Weight Polyethylene	XL-UHMWPE	ISO 5834
BIOLOX [®] <i>delta</i> ceramic		ISO 6474-2

Symbols

48 = 48mm spherical diameter, M7 = 7mm threaded hole.

Intended use

The ACE[®] Acetabular Cup System is a multi-component acetabular cup system designed for use without cement and in conjunction with a corresponding JRI femoral stem and femoral head, for primary or revision total hip arthroplasty.

Target Population

ACE[®] Acetabular Cup System is designed for use in skeletally mature individuals where there is sufficient bonestock to support implantation and fixation of the device.

Description

ACE[®] Acetabular Cup System comprises the following components:

ACE[®] H-A.C. Cup Shell is a true hemispherical acetabular shell made from Titanium alloy Ti₆Al₄V coated with Hydroxyapatite Ceramic, CA₅OH(PO₄)₃. The shells are designed to be implanted into the acetabulum without cement, fixation being achieved through bone ingrowth and union of the coating and the host bone. Various sizes of acetabular shells are available to accommodate anatomical variations of the acetabulum. Smaller sized implants are intended for patients with small bone and normally slight weight and could be inappropriate for other patients.

Screw holes permit the use of low profile titanium alloy bone screws for additional fixation and security, particularly in those cases where the acetabular bonestock is deficient. Shells are designed to minimize effects of osteolysis by utilizing screw holes that can be plugged. All cups are supplied with coated plugs fitted; these should be removed if a screw is to be used.

ACE[®] Low Profile Bone Screws are available for use to enhance fixation of the ACE[®] H-A.C. Cup Shell when required. They come in a variety of lengths and require use of a pre-drill manufactured of a non-implant grade stainless steel. This drill is supplied sterile and is intended for single use only.

Ensure that the heads of cancellous bone screws used with Acetabular Cups DO NOT PROTRUDE above the bore tapered surface before inserting the liner.

ACE[®] H-A.C. Cup Zero BIOLOX[®] *delta* Ceramic Liners are designed to sit within the acetabular cup and are fixed into the cup with an 18° taper lock. The liners are available in various sizes to accommodate different shell sizes and femoral head sizes. The ceramic liners should only be used with ceramic femoral heads; never use a metal femoral head with a ceramic liner. Special Instructions for use are required for ceramic heads and liners – refer to JRI Instructions For Use 155-020.

ACE[®] H-A.C. Cup CLP75 xl-UHMWPE Liners & Hooded Liners are made from highly cross-linked UHMWPE for further reduction in wear. They are designed to sit within the acetabular cup and are positioned by aligning the tabs on the rim of the liner with the 12 castellations on the rim of the shell. The poly liners have a secondary locking feature in the form of a ring which clips into the internal groove of the shell to prevent subluxation of the liner. A range of 10° and 20° hooded liners are available where there are concerns of dislocation. It should be noted that use of these liners slightly decreases the range of motion. The liners are available in various sizes to accommodate different shell sizes and femoral head sizes. XI-UHMWPE liners can be used with either a ceramic or metal femoral heads.

ACE[®] H-A.C. Cup Dual Mobility Insert comprises a CoCr sleeve and an ACE Dual Mobility xl-UHMWPE bearing to offer an uncemented Dual Mobility option, where there are concerns of dislocation. The liners are available in various sizes to accommodate different shell sizes and femoral head sizes. This dual mobility option must only be used in conjunction with either a cobalt chrome or a ceramic head.

Instructions for Use are required for ceramic heads and liners: Refer to JRI Instructions for Use 155-020

NOTE: the ACE[®] Acetabular Cup System components are NOT interchangeable with components of the other Furlong[®] H-A.C. Acetabular Cups, namely the CSF and CSF Plus.

Compatibility

The **ACE[®] Acetabular Cup System** is only compatible to use with JRI Orthopaedics Ltd. femoral stems and femoral heads. These include:

- **Furlong[®] H-A.C. Stems and Revision Stems**
- **Furlong[®] Modular Cemented Stem and Revision Stem**
- **Furlong[®] Evolution Stem**
- **Securus Revision Stem System**
- **AVANTEON Cemented Stem**

The aforementioned stems have their own instructions for use, which can be accessed via the JRI Orthopaedics Ltd. website or requested via the local sales representative.

Indications

The **ACE[®] Acetabular Cup System** is indicated for use in the following conditions:

- Severely painful and/or disabled joint from osteoarthritis, traumatic arthritis, rheumatoid arthritis, or congenital hip dysplasia
- Avascular necrosis of the femoral head
- Acute traumatic fracture of the femoral head or neck
- Failed previous hip surgery including joint reconstruction, internal fixation, arthrodesis, hemi-arthroplasty, surface replacement arthroplasty, or total hip replacement
- High risk of dislocation (**ACE[®] H-A.C. Cup Dual Mobility**)

Contra-indications

The device should NOT be implanted where there is active infection, insufficient bonestock to support the prosthesis or provide adequate fixation.

Further contra-indications may include the following conditions:

- Severe muscle, nerve or vascular diseases that endanger the extremity in question
- Severe deformations, tumours
- Severe osteoporosis or deficient bone substance that may endanger stable seating of the prosthesis
- Systematic and metabolic disorders

- Weakened or compromised immune system (HIV, tumours, infections)

Precautions

The following conditions require caution and due consideration during pre-operative planning by the surgeon:

- Obese or severely overweight patients
- Excessive loading through arduous activity
- Lack of mental faculties to understand post-operative recuperative regime
- Alcohol dependency or drug abuse
- History of falls or disabilities

In patients' with a high body mass index, where delayed surgery is feasible, it is advisable that a programme of weight reduction is undertaken prior to any joint replacement surgery.

Where the use of the dual mobility system is planned, patients should be assessed for **nickel sensitivity**, as the raw material used for the manufacture of some of the components may leach nickel in concentrations sufficient to elicit an allergic response in sensitive people.

Warning

This device should only be implanted by operating surgeons who are familiar with the general problems of prosthetic surgery and who are sufficiently trained to the product-specific operative technique. The surgeon is responsible for ensuring the surgery is carried out properly and in line with the instructions provided in the operative technique.

As a manufacturer, JRI Orthopaedics Ltd. Is not responsible for any complications arising from incorrect diagnosis, incorrect choice of implant, incorrect operative technique, treatment methods limitations or inadequate asepsis.

Pre-operative

Pre-operative planning allows the surgeon to assess for implant size and restoration of biomechanics. X-ray templates should be used to estimate implant sizes, placement and joint alignment. Failure to carry out proper planning may lead to incorrect choice of implant type/size. Ensure all implant sizes and required instrumentation are available prior to surgery. Consult operative technique and training materials provided by JRI Orthopaedics Ltd. before use.

Functionality of surgical instruments should be checked. Use of damaged instruments may lead to early failure of implant.

Transient bacteraemia can occur after surgical procedures. To prevent late infection at the implant site, many orthopaedic surgeons advise the use of antibiotic prophylaxis before and after such procedures for their patients.

Intra-operative

Always use a trial for any test fit and to check the range of motion. Failure to use the optimum size of implant of ensuring sufficient fixation may result in early revision surgery.

ACE[®] H-A.C. Acetabular Cups

Seat the cup at 45° inclination and 20° anteversion for proper positioning to decrease the chance of dislocation.

Ensuring that the cup is seated without the use of screws requires a greater level of surgical experience and it is recommended that experience is first gained using screws. The optimum number of bone screws should be used with the ACE Cups, so as to provide adequate stability, whilst minimising any fretting effects or the potential of screw failure.

It is advisable to pre-drill holes for screws using the provided Single Use drills and drill guide. Care should be taken not to over tighten bone screws. Failure to use the drill guide correctly may result in broken drills. Ensure the

appropriate selection of bone screw length and location to avoid damage to underlying soft tissue. Following screw insertion ensure that the screw heads are properly seated below the inner surface of the cup.

Clean surgical debris from the interior of the shell prior to seating the liner to prevent inadequate liner/shell interlock. Failure to properly seat the liner into the cup can lead to liner separation. Always ensure proper alignment and seating of the acetabular liner before impaction to prevent damage.

Avoid damage to the HA coating which could lead to debris particles. Implants MUST NOT be re-used under any circumstances as the functionality, integrity and/or sterility of that device may have been adversely affected and therefore cannot be guaranteed. Implant components from JRI Orthopaedics Ltd. should not be used together with those of another manufacturer, since compatibility of mating parts cannot be assured.

ACE® H-A.C. Cup Dual Mobility – When dual mobility components are to be used in conjunction with the **ACE® H-A.C. Acetabular Cups**, an intermediary cobalt chrome sleeve must first be fitted in the cup, to provide the necessary secondary bearing with the liner bearing. The bearing surface on the cobalt chrome **ACE® Dual Mobility** sleeve, when used in conjunction with the cementless acetabular cup, must be free of any debris, so as not to restrict movement between the liners and the bearing surfaces.

In addition, the Dual Mobility head assembling tool must be used when assembling both the trial and definitive implant dual mobility bearings and corresponding femoral heads. This is achieved by press-fitting the femoral head into the dual mobility bearing.

Warning

Care must be taken when determining and selecting the proper length of screws. Complications including internal bleeding and damage to vital organs have been reported as a result of trans-acetabular pelvic penetration. The risk of injury to vascular or neuralgic structures may be reduced by placing screws in the posterior inferior and posterior superior quadrants of the acetabulum.

Post-Operative

Physicians should ensure patients are aware of implant loading limitations and ensure consistent post-operative care, from a suitably qualified professional, is made available. The incidence and severity of complications are usually greater in surgical revisions than primary operations. Early detection of an impending complication allows for timely and effective countermeasures.

Post-operative care should incorporate recognised procedures and take into account information from the operative technique and documented according to internal hospital procedures. Failure to do this may result in malalignment, delayed wound/bone healing, implant failure, infection or impaired joint function.

Side Effects

As with all major surgical procedures, side effects and adverse events can occur. Some of the more common complications include:

- Problems resulting from anaesthesia and patient positioning (pain, nausea etc.)
- Infection
- Aseptic loosening of implant
- Dislocation, subluxation, insufficient range of movement, leg length discrepancy
- Damage to soft tissue
- Fracture of implant or bone
- Allergy/hypersensitivity reactions

- Cardiovascular/pulmonary embolism, respiratory infection, venous thrombosis, neuronal dysfunction, haematoma or delayed wound healing

Revision

CSF Plus Multi-hole Revision Cups and larger **ACE® Acetabular Cup** are available for use where bonestock may be compromised and feature larger sizes and a greater range of holes to provide the optimum position of the screw in the best bonestock. Fixation and press-fit of the cup should be considered to provide adequate fixation in the poorer quality bone.

Ensure all fragments of the primary prosthesis and any bone cement (if applicable) are removed. Clean and prepare the area for implantation in line with current Instructions for Use and Operative Technique of selected revision implant, prior to insertion of device.

If an UHMWPE liner is to be replaced, the acetabular cup can remain *in situ*, provided that there is still adequate fixation to the acetabulum

Special instructions apply when ceramic components are revised. In revision of any ceramic component, always use a ceramic formal head. If a ceramic liner is to be replaced, the acetabular cup will also need to be removed, as ceramic fragments may have damaged the inner surface/taper of the cup. New ceramic liners should **NOT** be fitted into a used cup taper. **Never** use a metal femoral head following the revision of any ceramic component. If revision is required in case of a fracture of a ceramic component, make sure all ceramic particles have been removed. Refer to JRI Instructions for Use 155-020 (ceramic components) prior to revision surgery for more details.

Storage & Handling

Store implants in their original protective packaging in a clean and dry atmosphere.

Do not use if the packaging is damaged. Do not use this product after the expiry date (year-month) shown on the product packaging.

Avoid removing from packaging until immediately prior to use. Inspect device prior to use. Visibly damaged, scratched, improperly handled implants and implants that have already been used must not be implanted under any circumstances as the functionality, integrity and/or sterility of that device may have been adversely affected and therefore cannot be guaranteed. Exposed articular surfaces must be neither marked nor come into contact with metallic or hard objects. Touching the articular surfaces, the Hydroxyapatite coating or the tapered interface for the liner must be avoided.

MRI Safety Information

Non-clinical testing has demonstrated that ACE® Acetabular Cup System are MRI Conditional. A patient with this device can be safely scanned in an MR system meeting the following conditions:

- Static magnetic field of 3.0 T
- Maximum spatial field gradient of 720 gauss/cm
- Maximum whole body averaged specific absorption rate (SAR) of 3.0W/kg

Under the scan conditions defined above, the **ACE® Acetabular Cup System** are expected to produce a maximum temperature rise of 1.7°C after 15 minutes of continuous scanning.

The image quality may be compromised if the area of interest is in the exact same area, or relatively close to, the position of the medical device. Therefore, optimization of MR imaging parameters to compensate for the presence of this device may be necessary:

Pulse Sequence	T1-SE	T1-SE	GRE	GRE
Signal Void Size	42,763-mm ²	17,875-mm ²	56,249-mm ²	28,917-mm ²
Plane Orientation	Parallel	Perpendicular	Parallel	Perpendicular

Additional information regarding MRI safety is available upon request at JRI Orthopaedics Ltd.

Functional Device Lifespan

The functional lifespan of an implant may be impacted by surgeon and selected operative technique, patient physiology and activity levels. Whilst it is normally expected that the lifetime of this implant will exceed a minimum of 10 years, it will be subject to wear and tear through normal everyday use.

Products intended for single-use must not be re-used

Re-use or reprocessing (e.g. cleaning and re-sterilisation) may compromise the structural integrity of the device and / or lead to device failure which may result in patient injury or death.

Reuse or reprocessing of single-use devices may create a risk of contamination (e.g. transmission of infectious material) which may result in injury or death.

Further information

For product specific operative technique training or further information, please contact your JRI Orthopaedics Ltd. Sales Representative or JRI Orthopaedics Ltd. directly.

JRI ORTHOPAEDICS LTD

18 Churchill Way,
35A Business Park,
Chapelton,
Sheffield,
S35 2PY, UK

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Ytterligare information avseende MRT-säkerhet finns tillgänglig på begäran hos JRI Orthopaedics Ltd.